


## STERIS VHP Integrated Systems

Discover effortless high-capacity, facility-wide decontamination

## WHAT IS VAPORIZED HYDROGEN PEROXIDE?



Broad Spectrum Anti-Microbial


Wide Range of Efficacy

High Level of Material Compatibility


Green Process

## WHAT IS AN INTEGRATED SYSTEM?

Integrated VHP Systems are installed into buildings so that the VHP generator is in direct communication with building management systems and air handling systems. VHP flows into rooms, isolators or other small enclosures using existing HVAC duct work or dedicated VHP piping.




**STERIS INTEGRATED VHP SYSTEMS ARE  
CONTROLLED BY BUILDINGS' AUTOMATION SYSTEM**





Our VHP Systems can be integrated into new construction or existing facilities.

## POWER AND VERSATILITY

Multiple areas may be decontaminated in sequence.



Preprogrammed cycles ensure repeatability of the calibrated cycles for each unique area, making operation simple and reliable.



At The  
Push of a  
Button

A Complete  
6 Log  
Cycle is Run

## — CHOOSE AN INTEGRATED SYSTEM FOR YOUR FACILITY SIZE —



### VHP M1000-T4

High capacity for facility wide decontamination



### VHP M100 & M100X

Unit for Isolators & Mid-Sized Enclosures



### VHP M10

Unit for Isolators & Small Enclosures

## VHP TURNKEY SERVICES

System Design | Integration Support | Cycle Development | Maintenance

## VHP CONSUMABLES



### Biological Indicators



### Vaprox® Hydrogen Peroxide Sterilant



### Process Indicators

Contact your STERIS Representative for additional information or to request a quote!